



**Easy application** for wide patient range.

# EMBOLECTOMY CATHETER



**Bioptimal™**

# EMBOLECTOMY CATHETER

Easy application  
for wide patient range.



- Bioptimal Embolectomy Catheter offer simple, rapid removal of emboli from the arterial system.
- Proprietary “low-profile” winding technique for balloon attachment insures balloon symmetry, providing uniform contact with vessel walls.
- Soft, rounded tip promotes easy arterial insertion for reduced trauma, arteriopuncture and plaque avulsion.
- Flexible polyvinylchloride body with “workhorse” latex balloon.
- Winged hub for easy catheter manipulation during insertion and attachment of balloon inflation syringe.
- Available in 4 lengths and 6 diameters, ranging from 2 French to 7 French.
- Size and balloon inflation capacity printed on proximal end for convenient and rapid identification.

## Specifications

Product Code	Catheter Size	Length (cm)	Colour Code	Maximum Liquid Volume (cc)	Maximum CO <sub>2</sub> Volume (cc)	Inflated Balloon Diameter (mm)	Maximum Pull Force Inflated Balloon (lb)
EL 3260	2F	60	Purple	0.08	0.2	4	0.5
EL 3280	2F	80	Purple	0.08	0.2	4	0.5
EL 3340	3F	40	Green	0.20	0.6	6	0.7
EL 3380	3F	80	Green	0.20	0.6	6	0.7
EL 1440	4F	40	Red	0.75	1.7	9	1.5
EL 1480	4F	80	Red	0.75	1.7	9	1.5
EL 14100	4F	100	Red	0.75	1.7	9	1.5
EL 1540	5F	40	White	1.50	3.0	11	2.0
EL 1580	5F	80	White	1.50	3.0	11	2.0
EL 1640	6F	40	Blue	2.00	4.5	13	2.5
EL 1680	6F	80	Blue	2.00	4.5	13	2.5
EL 1740	7F	40	Yellow	2.50	5.0	14	3.5
EL 1780	7F	80	Yellow	2.50	5.0	14	3.5

\*\*\*Please contact your Bioptimal representative for assistance.

### BIOPTIMAL INTERNATIONAL PTE. LTD.

36 Jalan Tukang  
Singapore 619266  
Tel: +65 6213 5777  
Fax: +65 6213 5737  
Email: sales@bioptimalg.com



[www.bioptimalg.com](http://www.bioptimalg.com)